

1 SECOND MEASUREMENT

Medical infrared thermometer XL-F02



1s



Precision probe



1 button measurement



3 colors backlight



32 sets of memories



Multi purpose



Infrared measurement

Product appearance



Ergonomic design




Accurate measurement

Sensor with exquisite craft
Temperature data collection



Precise probe
Smart calibration
Sensor with exquisite craft




Smart three colors back light reminder

According to the color of LCD screen, people can know the body temperature easily.



LCD digital screen

Digital LCD large screen, understand fully at one glance easier reading.



Product specification

Product name:	Medical infrared thermometer
Product model:	XL-F02
Product type:	Non-contact infrared thermometer
Measuring distance:	5cm-10cm
Measurement time:	≤ 2s
Automatic shutdown:	≤ 20s
Memory data:	32 sets
Measurement range:	32.0°C-42.9°C (89.6°F-109.22°F) (Body temperature mode)
Weight and size:	155 x 115 x 45 mm